


MALTOSE PH JP, CP

For Culture Media

What is MALTOSE PH ?

- Maltose is a reducing disaccharide consisting of two glucose molecules linked by an α -1,4 bond
- Maltose is manufactured from starch by enzyme technology
- Maltose PH is highly purified crystalline maltose monohydrate with low endotoxin

Structure



Chemical Formula: $C_{12}H_{22}O_{11} \cdot H_2O$


Molecular weight: 360.31

CAS RN[®]: 6363-53-7

Effect of Maltose PH supplementation on antibody production by CHO cell

Methods

- CHO-K1 cells were inoculated at 0.3×10^6 cells/mL in duplicates in a DMEM/F12-based protein free chemically defined medium (PFCDM) supplemented with 2 g/L D-(+)-glucose, or glucose with an additional 10 g/L MALTOSE PH in single-use Erlenmeyer flasks.
- The cultures were incubated in a humidified incubator at 37°C, 8% CO₂ and a rotation speed of 110 rpm. Cell culture supernatants were collected daily on days 0-5 and 7, and monoclonal IgG antibody (anti-Her2) titers were determined by nephelometry using IMMAGE 800 (Beckman Coulter).



Results

- Two g/L of glucose was chosen as the base glucose concentration because this will allow for a premature but controlled cell growth constraint due to glucose depletion on Day 4.
- IgG titer reached a maximum of 340 mg/L on Day 5 in the culture containing 2 g/L glucose as a carbon source.
- IgG production continued in the MALTOSE PH supplemented cultures Days 5 through 7 to reach 600 mg/L.

Manufacturer : **HAYASHIBARA CO., LTD.**

CONTACT : **NAGASE & CO., LTD.**

Life & Healthcare Products Dept. Pharma-Medical Div.

TEL: +81-3 (3665) 3333 (TOKYO JAPAN)

TEL: +81-6 (6535) 2327 (OSAKA JAPAN)

E-mail: dnfct@ex.nagase.co.jp

The information provided herein is intended only for reference purposes. It is the customer's responsibility to determine that the ingredient meets all legal requirements in the country where it is used, and that it does not infringe on any third party patents.


Unauthorized reproduction of this brochure is prohibited.

For Culture Media

Effect of Maltose supplementation on cell viability and viable density

Methods

- The CHO-K1 cells were inoculated and cultured as described in the previous experiment.
- Viable cell density (VCD) and culture viability were analyzed by Vi-Cell XR Cell Viability Analyzer (Beckman Coulter, Brea, CA) according to manufacturer's instructions.



Results

- Growth of CHO-K1 cells reached a maximum viable cell density (VCD) of 8.9×10^6 cells/mL on Day 5 when cultured using 2 g/L glucose.
- MALTOSE PH supplemented cultures attained a higher maximum VCD of 10.2×10^6 cells/mL and longer culture viability compared to the culture containing only 2 g/L glucose.

All the data were provided by Dr. Say Kong Ng of Bioprocessing Technology Institute (BTI)

Summary

- MALTOSE PH can be used by the CHO-K1 cells as a carbon source to maintain culture viability and IgG production upon glucose depletion.
- MALTOSE PH supplementation to the CHO-K1 cell culture in addition to glucose increased the IgG titer compared to the culture with only glucose as a carbon source.
- MALTOSE PH supplementation along with glucose improved the growth of CHO-K1 cell compared to the culture without MALTOSE PH supplementation.
- It demonstrates that MALTOSE PH can effectively improve the productivity of antibody when used to supplement glucose as a carbon source.

Reference

Application of maltose as energy source in protein-free CHO-K1 culture to improve the production of recombinant monoclonal antibody. Leong DSZ *et al.* Sci Rep. 2018 Mar 6;8(1):4037.

Product Information

HIGH PURITY MALTOSE

Product Name	Purity/Other	Packaging	Regulatory Approval	Others
MALTOSE PH	Not less than 98.0% /Low Endotoxin	25kg PE bag in a carton box	JP	<ul style="list-style-type: none"> • US Type II DMF • Halal

Manufacturer : **HAYASHIBARA CO.,LTD.**

CONTACT : **NAGASE & CO.,LTD.**

Life & Healthcare Products Dept. Pharma-Medical Div.

TEL: +81-3 (3665) 3333 (TOKYO JAPAN)

TEL: +81-6 (6535) 2327 (OSAKA JAPAN)

E-mail: dnfct@ex.nagase.co.jp

The information provided herein is intended only for reference purposes. It is the customer's responsibility to determine that the ingredient meets all legal requirements in the country where it is used, and that it does not infringe on any third party patents.

Unauthorized reproduction of this brochure is prohibited.