

PEARLITOL® 100 SD - MANNITOL

Definition

Product Identifier

Product name: PEARLITOL® 100 SD - MANNITOL

D-MANNITOL

INCI: MANNITOL

Specifications

A) CHARACTERS

APPEARANCE White, crystalline powder

B) IDENTIFICATION

IDENTIFICATION TEST-A	IP	COMPLIES
IDENTIFICATION TEST-B(*)	IP	COMPLIES
IDENTIFICATION TEST-C(*)	IP	COMPLIES
IDENTIFICATION TEST-D	IP	COMPLIES

C) TESTS

D-MANNITOL (on DS)	IP	97.0 to 102.0 %
APPEARANCE IN SOLUTION	IP	clear, colourless
ACIDITY (NaOH 0.01M)(*)	IP	0.2 mL max.
SPECIFIC OPTICAL ROTATION(*)	IP	+23.0 to +25.0°
ARSENIC(**)	IP	2 mg/kg max.
CHLORIDES	IP	50 mg/kg max.
SULPHATES	IP	100 mg/kg max.
REDUCING SUGARS(Na ₂ SO ₃)	IP	12.8 mL min.
SORBITOL(*)	IP	COMPLIES
SULPHATED ASH	IP	0.1 % max.
LOSS ON DRYING	IP	0.5 % max.
PARTICLE SIZE (laser)		
- PART. > 150 microns (BECKMAN)		35 % max.
- PART. > 140 microns (SYMPATEC)		35 % max.
- PART. > 75 microns (BECKMAN)		70 % min.
- PART. > 30 microns (SYMPATEC)		70 % min.

PRODUCT SPECIFICATIONS SHEET

PEARLITOL® 100 SD - MANNITOL

MICROBIOLOGICAL VALUES:

- TOTAL AEROBIC MICROBIAL COUNT	1000 CFU/g max.
- TOTAL YEASTS AND MOULDS COUNT	100 CFU/g max.
- ESCHERICHIA COLI(**)	Not detected in 1g
- SALMONELLA(**)	Not detected in 10g

- * Compliance data - Tests not performed
- ** Monitoring plan

Comments

- Not intended for use in manufacture of parenteral dosage forms.

- NOTE FOR LASER PARTICLE SIZE DESCRIPTION

As part of continuous improvement, ROQUETTE decided to replace its current laser particle size device BECKMAN brand by the new technology SYMPATEC brand. To support this evolution, the specification has slightly evolved for particle size description, the acceptance criteria remain the same and the sieve size shows the both technology brand name and equivalency. The new expression with SYMPATEC technology is shown in this product specification sheet. ROQUETTE is currently in a transition period during which batches can be tested by BECKMAN device or SYMPATEC device.

Conformity

Conforms to the requirements of the current monograph

- Indian Pharmacopoeia (IP) MANNITOL

Please contact us for any statement regarding compliance to the General Chapters (elemental impurities, residual solvents, organic volatile impurities, metal catalyst, metal reagent).

Methods used by Roquette may be the Pharmacopoeia methods or alternative validated methods which have been compared to the Pharmacopoeia methods.

Storage

Expiry date Manufacturing date + 5 years, in its unopened packaging.

- The product durability may vary according to packaging type and manufacturing plant. Proper information is shown on labelling and CoA.
- We recommend to preserve the product in its unopened original packaging, preferably protected from wide variations of temperature and humidity.
- Upon opening, use the product as quickly as possible to prevent moisture regain.

Disclaimer



PEARLITOL® 100 SD - MANNITOL

The information provided in this Product Specification Sheet relates only to the specific product designated and may not be applicable when such product is used in combination with other materials or in any process.

All information and instructions provided in this Product Specification Sheet are based on the current state of our knowledge at the latest revision date indicated. It is the responsibility of the user to be aware of and to follow the regulations applying to our product for its possession, handling and use.

Notes : All the dates are formatted like YYYY/MM/DD.